



Knox® MedVault®



HIGH SECURITY NARCOTICS LOCKER

Controlled substance security is an ever-increasing concern for both fire departments and independent EMS providers. Drugs such as Morphine and Valium require top security, yet emergency personnel need immediate access during an emergency.

Paramedics can now significantly reduce the risk of carrying controlled substances by installing the new Knox MedVault® Narcotics Locker.

The vault is made of heavy aluminum and secured with a high security electronic keypad. Personnel gain access by entering their unique PIN code.

All activity at the unit is collected in an audit trail. This audit trail is collected via a USB port and provides a time and date stamp, including a user ID for each activity. This provides the department accountability regarding access to controlled substances.



Surface Mount Model

- Unit opens with PIN Code authorization
- Stores 1 to 5,000 unique PIN codes and 25,000 events
- Connects to PC USB port to collect release history (audit trail events)
- Optional Remote Administration
- Operates with a single 12 volt power source
- Light weight construction designed for use in mobile emergency/rescue vehicles



(Bag not included)

Heavy-duty bolt action provides the security you've come to expect from Knox




Recess Mount Model
MedVault® is available with or without recess mount flange

Standard Hirose USB Data Port
For audit trail collection

Power Cable

Optional Rear USB Data Cable
For Remote Administration

(cable preinstalled on models 5503 and 5504)



Size:	9-1/4" H x 12" W x 9" D Recess mount: Flange protrudes 7/8" on all sides Usable interior space: 7-1/2" H x 10-1/2" W x 6-1/2" D
Weight:	18 lbs
Finish:	Powder Coating
Color:	Grey
Material:	Heavy Aluminum, all welded construction
Features:	<ul style="list-style-type: none"> • Great strength with light weight construction for mobile use consideration • PIN code access (one or two PIN Codes for release) • Manual Key override in case unit loses power • Audit trail available via USB port



Knox Company
1601 W. Deer Valley Road
Phoenix, AZ 85027

Tel 800-552-5669
Fax 623-687-2299

www.knoxbox.com
E-mail: info@knoxbox.com





HIGH SECURITY NARCOTICS LOCKER

MedVault® Software

Knox MedVault software allows a department to manage all PINs needed to access the vault. Additionally, the audit trail history is accessed using the application software. All MedVault files are encrypted and require a password to access.



Each MedVault unit comes standard with a Virtual Interlock Validator (VIV). This electronic security feature requires the door to be closed prior to locking. For added security the encrypted serial number of the VIV is unique to each unit.

MedVault software is written for emergency responders and allows a department to program their units for varying degrees of security to meet the unique needs of their department.

MedVault® Programming Options

- **An individual with a valid PIN can open the vault; PIN is active 24/7**
- **Two individuals must supply valid individual PINs to open the vault; PINs are active 24/7**
- **An Individual with a valid PIN can open the vault; PIN is active only when shift is on-duty**
- **Two individuals must supply valid individual PINs to open the vault; PINs are active only when shift is on-duty**
- **Other security scenarios available; contact your local Knox account manager for details**

Optional Remote Administration Available

- **Program and collect data remotely**
- **Use city-wide or departmental WiFi network for communication with MedVault units**
- **Perform administrative duties from a central location:**
 - **Firmware Updates**
 - **Audit trail collection**
 - **PIN Code Changes, Additions and Deletions**

Coming Soon – Mini MedVault®



Standard MedVault

Mini MedVault

The Mini MedVault provides the same access control as our standard MedVault. Controlled substances are secured until authorized personnel access them with their unique PIN code. An audit trail showing the time, date and PIN code is recorded each time the vault is opened. This audit trail is collected via a USB port or from the convenience of an office with optional Remote Administration. This unit is specifically designed and constructed for mounting in a rescue vehicle or aircraft.

Mini MedVault

Dimensions: Surface: 6-3/4”H x 9”W x 5-1/4”D
Weight: 12 lbs
Finish: Powder Coating
Color: Grey

Check with your local account manager for availability and pricing